

# ***DuraTray II***

## **Custom Tray Resin**

**DuraTray II is a truly amazing  
and innovative custom tray resin  
that can help you cut your  
working time in half**

This system has been specially designed to produce a non-sticky, incredibly adaptable smooth dough in about 30 seconds. You will notice a substantial reduction in fabrication time since the dough is workable right up to the cure point in about 6 minutes.



*Since 1894: Over 100 years  
of Reliable Service to the Dental Community*

This system is designed to produce a non-sticky, incredibly adaptable dough in about 30 seconds. It minimizes fabrication work, is workable right up to the cure point, and is constantly self-adhesive. A handle can be added at any point without having to add monomer or mix additional resin.

A 1 lb. package contains a 16 oz. powder, 4 oz. liquid, a measuring vial, mixing cup, spoon, and dropper.

#### **INSTRUCTIONS:**

1. Pour liquid into orange vial. Fill to top of vial and pour into mixing cup.
2. Pour powder into blue spoon, slightly underfilling. The powder should be just below the top of the spoon. Pour into cup.
3. Mix thoroughly for 30 seconds.
4. **DuraTray II** is now ready to adapt.



**Dental Mfg. Co.**

5805 West 117th Place  
Worth, IL 60482

(708) 597-6694

FAX (708) 597-7560