



Reliance Dental Manufacturing, LLC

Why you should be a **Dura** Dentist



Reliance Dental Manufacturing, LLC, provides only the highest quality products, affordable prices, and personal attention to every customer.

Our **Dura** line of dental acrylics features outstanding accuracy, reliability, and ease of use. We take pride in knowing our products are designed to provide both **comfort** and **value** — ensuring a better experience for both patient and dentist.



Dura ... the first name in reliable dentistry.

DuraLay — Our signature inlay pattern resin sets in 5 minutes, forming a dense, accurate, non-shrinking pattern that can be trimmed, handled and invested without fear of breakage, chipping, or distortion.

DuraLay Temporary Crown & Bridge — A self-curing tooth shaded acrylic. Temporary Crown & Bridge provides very little shrinkage, is very hard, and easily polished and trimmed. Sets up in 4 to 5 minutes.




DuraSeal — A self-curing tooth shaded acrylic used as a temporary filling to seal cavities, temporary inlays or crowns and bridges.

DuraBase Soft Cushion Rebase — Self-curing, soft-cushion rebase material for relining dentures. Will last 9-12 months in the mouth. Provides relief and comfort for even the most sensitive patients.

Dura Rely-A-Soft — A long-term soft denture liner for chairside or laboratory use that provides added comfort and will last up to 2 years in the mouth. Available in pink or clear.

DuraLiner II — A permanent, hard, self-curing chairside rebase material for relining dentures. Low exothermic reaction; will not sting or burn the patient's tissue. Available in pink, clear and nylon fibers.

DuraConditioner — A tissue conditioner and functional impression material that treats and tones irritated tissue. Generates a functional impression.



See our full line of dental acrylics and supplies at **RelianceDental.net**

Over 100 years of Reliable Service to the Dental Community

Reliance Dental Manufacturing has been one of America's most trusted manufacturers of acrylics and dental supplies for generations. Our products are used daily in countries across the globe providing the quality, value, and customer satisfaction dentists can rely on.



Reliance Dental Manufacturing, LLC

5805 W. 117th Place • Alsip, IL 60803
(708) 597-6694
RelianceDental.net